## **☑** TriDFusion (3DF) Image Viewer KEYBOARD SHORTCUTS

OPERATION	KEYBOARD SHORTCUT	ICON	OPERATION	KEYBOARD SHORTCUT	ICON
PRESET WINDOW/LEVEL			DRAW CONTOUR		
CT windows	F1 (1200, -500) Lung		Draw freehand	D	
	F2 (500 , 50 ) Soft		Draw polygone	P	
	F3 (500 , 200 ) Bone		Draw sphere	E	<b>(</b>
	F4 (240 , 40 ) Liver		Click-VOI	V	G
	F5 (80 , 40 ) Brain			2D BRUSH & TOOL SET	
	F6 (350 , 90 ) Head and Neck		2D brush	В	, <b>,</b>
	F7 (2000, -600) Enhanced Lung		Change 2D brush diameter	Ctrl + Left Mouse Button Pressed Up/Down	
	F8 (350 , 50 ) Mediastinum		Knife	К	•
	F9 is assign to 3 Windows\Level:		Interpolate	Х	
	1. First press (2000, 0) All		VOI CONTOUR REVIEW		
	2. Second press (2500, 415 ) Vertebra		Select for contour review	Ctrl + S (with hovering over a ROI)	
	3. Third press (1000, 350) Temporal Bone		Show/Hide contours	S	
NAVIGATION CONTROLS			Add a freehand to a VOI	tab	
Zoom	Middle Mouse Button Pressed Up/Down		Previous VOI	, (<)	
Activate Zoom	Z (right click zoom off)	1	Delete VOI	delete	
Pan	Right Mouse Button Pressed Left/Right/Up/Down		Next VOI	. (>)	
Activate Pan	N (right click pan off)		Increase VOI contours size	Ctrl + +	
Fusion Alpha On/Off	The 'F' key is assigned to 3 operations:		Decrase VOI contours size	Ctrl + -	
	1. First press shows the fused series		ROI UTIL		
	2. Second press shows the reference series		Copy contour	Ctrl + C (with mouse over the ROI to copy)	
	3. Third press return to fusion		Paste contour	Ctrl + V	
Scroll	Shift + Left Mouse Button Pressed Up/Down		VOI = Volume-of-interest		
	Mouse Scroll Wheel		ROI = Region-of-interest		
Zoom in	+	⊕ <b>(</b>	DISPLAY OPTIONS		
Zoom out	-	Q	Show/Hide crosshairs	С	
Next slice	<b>↑</b>		Change crosshair type	Shift + C	
Previous slice	↓		Show/Hide overlay text	0	
Rotate MIP left	<b>←</b>	0	Invert color table	U	
Rotate MIP right	→	0	Bilinear Interpolation On/Off	I	
Full Screen Enter/Exit	Press Ctrl + F (while hovering over a view)	22			
	ORIENTATION & FLIP				
Flip Anterior/Posterior	A			CANCEL & RESET	
Flip Left/Right	L		Cancel Operation	Escape	
Flip Slices Order	R	7	Reset View	Space bar	7