

TriDFusion (3DF) Image Viewer KEYBOARD SHORTCUTS

OPERATION	KEYBOARD SHORTCUT
CT windows	F1 (1200, -500) Lung F2 (500 , 50) Soft F3 (500 , 200) Bone F4 (240 , 40) Liver F5 (80 , 40) Brain F6 (350 , 90) Head and Neck F7 (2000, -600) Enhanced Lung F8 (350 , 50) Mediastinum F9 is assign to 3 Windows\Level: First press (2000, 0) All Second press (2500, 415) Vertebra Third press (1000 , 350) Temporal Bone
Show/Hide crosshairs	c
Change crosshair type	shift+c
Show/Hide overlay text	o
Invert color table	i
Flip Anterior/Posterior	a
Flip Left/Right	l
Flip Slices Order	r
Show/Hide Contours	s
Draw freehand	d
Draw polygone	p
Draw sphere	e
Click-VOI	v
2D brush	b
Knife	k
Interpolate	x
Add a freehand to a VOI	tab
Change 2D brush diameter	Ctrl + Left Mouse Button Pressed Up/Down
Zoom	Middle Mouse Button Pressed Up/Down
Activate Zoom	z, right click zoom off
Pan	Right Mouse Button Pressed Left/Right/Up/Down
Activate Pan	n, right click pan off
Fusion Alpha On/Off	The 'f' key is assigned to 3 operations: First press shows the fused series. Second press shows the reference series. Third press return to fusion .
Scroll	Shift + Left Mouse Button Pressed Up/Down
Zoom in	+
Zoom out	-
Next slice	↑
Previous slice	↓
Rotate MIP left	←
Rotate MIP right	→
Cancel operation	Escape
Reset View	Space bar