TriDFusion (3DF) Image Viewer KEYBOARD SHORTCUTS

OPERATION	KEYBOARD SHORTCUT
CT windows	F1 (1200, -500) Lung
	F2 (500 , 50) Soft
	F3 (500 , 200) Bone
	F4 (240 , 40) Liver
	F5 (80 , 40) Brain
	F6 (350, 90) Head and Neck
	F7 (2000, -600) Enhanced Lung
	F8 (350, 50) Mediastinum
	F9 is assign to 3 Windows\Level:
	First press (2000, 0) All
	Second press (2500, 415) Vertebra
	Third press (1000 , 350) Temporal Bone
Show/Hide crosshairs	С
Change crosshair type	shift+c
Show/Hide overlay text	О
Invert color table	i
Flip Anterior/Posterior	a
Flip Left/Right	I
Flip Slices Order	r
Draw freehand	d
Draw polygone	р
Draw sphere	s
Click-VOI	v
2D brush	b
Knife	k
Add a freehand to a VOI	tab
Change 2D brush diameter	Ctrl + Left Mouse Button Pressed Up/Down
Zoom	Middle Mouse Button Pressed Up/Down
Activate Zoom	z, right click zoom off
Pan	Right Mouse Button Pressed Left/Right/Up/Down
Activate Pan	n, right click pan off
Fusion Alpha On/Off	The 'f' key is assigned to 3 operations:
	First press shows the fused series.
	Second press shows the reference series.
	Third press return to fusion .
Scroll	Shift + Left Mouse Button Pressed Up/Down
Zoom in	+
Zoom out	-
Next slice	<u>↑</u>
Previous slice	<u> </u>
Rotate MIP left	(
Rotate MIP right	→
Cancel operation	Escape
Reset View	Space bar