TriDFusion (3DF) Image Viewer KEYBOARD SHORTCUTS

OPERATION	KEYBOARD SHORTCUT	ICON	OPERATION	KEYBOARD SHORTCUT	ICON
PRESET WINDOW/LEVEL			DRAW CONTOUR		
CT windows	1 (1200, -500) Lung		Draw freehand	D	
	2 (500 , 50) Soft		Draw polygone	P	
	3 (500 , 200) Bone		Draw sphere	E	\oplus
	4 (240 , 40) Liver		Click-VOI	V	G
	5 (80 , 40) Brain			2D BRUSH & TOOL SET	
	6 (350 , 90) Head and Neck		2D brush	В	**
	7 (2000, -600) Enhanced Lung		Change 2D brush diameter	Ctrl + Left Mouse Button Pressed Up/Down	
	8 (350 , 50) Mediastinum		Knife	K	•
	9 is assign to 3 Windows\Level:		Interpolate	X	
	1. First press (2000, 0) All			VOI CONTOUR REVIEW	
	2. Second press (2500, 415) Vertebra		Select for contour review	Ctrl + S (with hovering over a ROI)	
	3. Third press (1000, 350) Temporal Bone		Show/Hide contours	S	
NAVIGATION CONTROLS			Add a freehand to a VOI	tab	
Zoom	Middle Mouse Button Pressed Up/Down		Previous VOI	, (<)	
Activate Zoom	Z (right click zoom off)		Delete VOI	delete	
Pan	Right Mouse Button Pressed Left/Right/Up/Down		Next VOI	. (>)	
Activate Pan	N (right click pan off)		Increase VOI contours size	Ctrl + +	
Fusion Alpha On/Off	The 'F' key is assigned to 2 operations:		Decrase VOI contours size	Ctrl + -	
	1. First press shows the fused series		ROI UTIL		
	2. Second press return to fusion		Copy contour	Ctrl + C (with mouse over the ROI to copy)	
			Paste contour	Ctrl + V	
	The 'G' key is assigned to 3 operations:		VOI = Volume-of-interest		
	1. First press shows the fused series		ROI = Region-of-interest		
	2. Second press shows the reference series			DISPLAY OPTIONS	
	3. Third press return to fusion		Show/Hide crosshairs	С	
Scroll	Shift + Left Mouse Button Pressed Up/Down		Change crosshair type	Shift + C	
	Mouse Scroll Wheel		Show/Hide overlay text	0	
Zoom in	+	⊕ (Invert color table	U	
Zoom out	-	Q	Bilinear Interpolation On/Off	I	
Next slice	↑				
Previous slice	↓				
Rotate MIP left	←	۵			
Rotate MIP right	→	0			
Full Screen Enter/Exit	Press Ctrl + F (while hovering over a view)	K X			
ORIENTATION & FLIP					
Flip Anterior/Posterior	A			CANCEL & RESET	
Flip Left/Right	L	_	Cancel Operation	Escape	
Flip Slices Order	R		Reset View	Space bar	